

Chennai Metco

for your metallurgy



Finest products for
Micro Structure Analysis



Consumables Catalogue

2016-17

www.chennaietco.com

SECTIONING





Choosing the Right Wheels for Right Application is important. Wrong choice of Wheels can sometimes lead to damaging your specimen.

Using Appropriate grade of wheels makes a big difference in the cutting performance.


The skill of manual cutting and appropriate cutting mode / feed combined with the right wheel is the recipe for success in abrasive cutting.

Breakage of wheels or damage to the cutting process can best be avoided by proper care.


ABRASIVE CUT OFF WHEELS

CUTTING WHEEL	APPLICATION	DIMENSIONS (mm)	CODE NO
GENERAL PURPOSE WHEELS 	Used on Ferrous Materials with Hardness upto 40 HRC	200 x 31.75 x 1	CW001R
		200 x 31.75 x 1.5	CW01R-1.5
		250 x 31.75 x 1.5	CW002
		300 x 31.75 x 2	CW003
		350 x 31.75 x 2.5	CW004
HARD STEEL WHEELS 	Used on Hard Ferrous Materials with Hardness 40 - 60 HRC	200 x 31.75 x 1	CW01S-1
		200 x 31.75 x 1.5	CW001S
		203 x 25.4 x 0.7	CW01S-0.7
		250 x 31.75 x 1.5	CW005
		300 x 31.75 x 2	CW006
NORTON WHEELS 	Coarse cutting wheel. Can be used after testing	300 x 25.4 x 2.5	CWN012
		350 x 25.4 x 3	CWN014
WHEELS FOR NON FERROUS & CAST IRON 	Non-ferrous Materials, Cast Iron	350 x 31.75 x 3	CWA008

BREAKAGE RESISTANT CUT OFF WHEELS

CUTTING WHEEL	APPLICATION	DIMENSIONS (mm)	CODE NO
PARTIALLY REINFORCED WHEELS 	Used on Ferrous Materials with Hardness up to 40HRC . Suited to Heavy duty components .	300 x 31.75 x 1.5	CW003PR
		400 x 31.75 x 3.2	CW009PR
		500 x 31.75 x 3.2	CW013P
		500 x 40 x 3.2	CW014
	Used on Hard Ferrous Materials with Hardness up to 60HRC . Suited to Heavy duty components	500 x 31.75 x 3.2	CW013S

HIGHER BREAKAGE RESISTANT CUT OFF WHEELS

CUTTING WHEEL	APPLICATION	DIMENSIONS (mm)	CODE NO
FULLY REINFORCED WHEELS 	Used on Ferrous Materials with Hardness up to 40HRC. More Breakage Resistant for Tough cutting applications	350 x 31.75 x 3	CWF008FR
		400 x 31.75 x 3.2	CW009FR
		500 x 31.75 x 3.2	CW015FR

ABRASIVE CUT OFF WHEELS

SUPREME CUTTING FINISH – USA GRADE





APPLICATION	DIMENSIONS (mm)	CODE NO
Used on <i>Abrasive cutting</i> Machines		
Tool Steels, Cast Iron (above 60HRC)	255 x 1.5 x 31.75	CW-031
Hardened Steels (above 50HRC)	255 x 1.5 x 31.75	CW-032
Medium Hard Steels and Super alloys (35-50 HRC)	255 x 1.5 x 31.75	CW-033
Medium Soft Metals (15-35 HRC)	255 x 1.5 x 31.75	CW-034
For Aluminium, Brass, Titanium and Zirconium	255 x 1.5 x 31.75	CW-036
Used on <i>Precision cutting</i> Saws		
For Titanium and other Alloys	152 x 1 x 16	647
General Ferrous	101 x 0.3 x 12.7	4F12
General Non Ferrous	101 x 0.3 x 12.7	4NF12
General Ferrous	127 x 0.5 x 12.7	5F20H
Soft to Medium hard Ferrous (below 45 HRC)	127 x 0.3 x 12.7	5F20S

**Other dimensions available*

**Rubber bonded wheels for high quality burn-free sectioning also available*

PRECISION CUT OFF WHEELS


(FOR SMALL COMPONENTS & THIN SLICING APPLICATIONS)- GERMAN GRADE

CUTTING WHEEL	APPLICATION	DIMENSIONS (mm)	CODE NO
 150mm wheels can be used in Precision Saws	Used on Mild Steel, Structural Steel, Free cutting steel, Spring Steel, Piano Wire, Perm alloy, Nichrome (Upto 40HRC)	150 x 0.5 x 25.4	PCW0001
		150 x 1.0 x 25.4	PCW0002
		205 x 0.8 x 31.75	PCW0003
		255 x 1.2 x 31.75	PCW0004
		305 x 1.2 x 31.75	PCW0005
 150mm wheels can be used in Precision Saws	Used on Hard Ferrous Alloys, Alloy Tool Steel, Bearing Steel, Stainless Steel, Heat resistant Steel, Invar Steel, Silicon Steel, Hardened Steels (Upto 60HRC)	150 x 0.5 x 25.4	PCW1001
		150 x 1.0 x 25.4	PCW1002
		205 x 0.7 x 31.75	PCW1003
		255 x 1.2 x 31.75	PCW1004
		305 x 1.2 x 31.75	PCW1005
		355 x 2.0 x 31.75	PCW1006
 150mm wheels can be used in Precision Saws	Used on Medium Soft Steels, Thin cut, Cast Iron, Cast Steel, KS Steel, Brass, Bronze, Aluminium alloys, Tungsten Molybdenum alloy, Pure Titanium(α alloy), Lead alloy, Glass (15-35HRC)	405 x 2.5 x 31.75	PCW1007
		150 x 0.5 x 25.4	PCW2001
		205 x 0.8 x 25.4	PCW2002
		230 x 1.2 x 31.75	PCW2003
		255 x 1.2 x 31.75	PCW2004


SPACER

PRODUCT	DESCRIPTION	SIZE	QTY	CODE NO
	Spacer to fit in 1" Spindle and wheel with 31.75 mm ID	25.4mm ID	10	CWS06

DIAMOND WAFERING BLADE

WAFERING BLADE	APPLICATION	DIMENSIONS (mm)	CODE NO
	Alumina ceramic, Silicon-Nitride ceramic, glass, copper, and iron, ruby, Quartz, Silica	100 x 1.27 x 0.3	T4DW
		150 x 31.75 x 1	CWD002
		200 x 31.75 x 1.2	CWD003
		250 x 31.75 x 1.4	CWD004

DIAMOND WAFERING BLADE USA GRADE

WAFERING BLADE	APPLICATION	DIMENSIONS (mm)	CODE NO
	Ceramics, Carbides, Boron Composites, Glass, Refractory (above 60HRC) Used in Precision Saws	75 x 12.7 x 0.15	M306L
		100 x 12.7 x 0.3	M412L
		125 x 12.7 x 0.4	M515L
		150 x 12.7x 0.5	M620L
		175 x 12.7x 0.6	M725L
	Polymers, PC Boards, Titanium, Thermal Spray Coating Used in Precision Saws	75 x 12.7 x 0.15	M306H
		100 x 12.7 x 0.3	M412H
		125 x 12.7 x 0.4	M515H
		150 x 12.7x 0.5	M620H
		175 x 12.7x 0.6	M725H
	Hard Ferrous materials	200 x 31.75 x 1.1	M845D
		250 x 31.75 x 1.3	M1050D
		300 x 31.75 x 1.6	M1262D
		355 x 31.75 x 1.8	M1470D

CBN WAFERING BLADE – USA GRADE

WAFERING BLADE	APPLICATION	DIMENSIONS (mm)	CODE NO
	High speed steels, Inconel, Stellite, Iron and Cobalt based alloys, Nickel based Super alloys, Lead base alloys, Titanium alloys and Stainless Steels. (Used on low speed /High speed)	75 x 12.7 x 0.15	M306C
		100 x 12.7 x 0.3	M412C
		125 x 12.7 x 0.4	M515C
		150 x 12.7 x 0.5	M620C
		175 x 12.7 x 0.6	M725C
	Used in Precision Saws High speed steels, Inconel, Stellite, Iron and Cobalt based alloys, Nickel based Super alloys, Lead base alloys, Titanium alloys and Stainless Steels. (Used on High speed)	200 x 31.75 x 1.1	M845C
		250 x 31.75 x 1.3	M1050C
		300 x 31.75 x 1.6	M1262C
		355 x 31.75 x 1.8	M1470C

DRESSING STICKS

PRODUCT	DESCRIPTION	SIZE	QTY	CODE NO
	Used on new and frequently used diamond wafering blades and grinding discs to maintained superior cutting performance	4" x 0.5"	1	DS-10

CUTTING FLUID


PRODUCT	QTY	WT Kg	CODE NO
COOLANT OIL	1 Ltr	1.2	CWCO07
Water soluble coolant, used in the ratio 1:30 (Oil:Demin Water) For abrasive cutting operations			
ISOMET OIL	1litre	1.2	DB101
Oil-based coolant which can efficiently reduce the Cycle time of Cutting. Used with precision cutting saws.			




MOUNTING

Mounting aids in easy handling of the sample and to preserve the edges of the sample.

Chennai Metco offers a wide range of mounting resins in hot and cold mounting systems that suits your application at its best.



COMPRESSION MOUNTING POWDERS


MOUNTING RESIN	COLORS	CODE NO	
PHENOLIC POWDER Used for regular analysis , with Fast cycle time	Black	MB101	
	Red	MB102	
	Green	MB103	
TRANSPARENT POWDER Suited to applications that require complete visibility of the sample, good flow characteristics	White	MT104	
TRANSPARENT GRANULES Suited to applications that require complete visibility of the sample, moderate flow characteristics	Transparent granules	MTG105	
CONDUCTIVE POWDER Applied across analysis that requires the a conductive sample , for SEM analysis	Brown	MB105	

PRODUCT NAME	DESCRIPTION	Qty	CODE NO
MOLD RELEASE SPRAY 	Sprayed on the walls of the mounting chamber for easy removal of the mold	200ml	MS001

For better edge retention choose premium quality products


COMPRESSION MOUNTING POWDERS

PREMIUM QUALITY

MOUNTING RESIN	COLORS	CODE NO	
PHENOLIC POWDER			
Used for Regular analysis with of Medium hard materials	Black	M126	
	Red	M127	
	Green	M128	
<i>Moderate Edge retention</i>			
EPOXY POWDER(COARSE GLASS FILLED)	Black	M129C	
Used for Regular analysis of Hard Materials			
<i>Good edge retention</i>			
EPOXY POWDER(FINE GLASS FILLED)	Black	M129	
Used for High end applications like coating thickness and case depth measurement, fastener profile analysis, etc			
Applied on Hard materials			
<i>Excellent edge retention</i>			
EPOXY POWDER(MINERAL FILLED)	Black	MP5000	
Used for High end applications like coating thickness and case depth measurement, fastener profile analysis, etc			
Applied on High Hardness materials			
<i>Excellent edge retention</i>			

COMPRESSION MOUNTING POWDERS

PREMIUM QUALITY

MOUNTING RESIN	COLORS	CODE NO	
Diallyl Phthalate (Glass Filled) Used for High end applications like coating thickness and case depth measurement, fastener profile analysis, etc Applied on Moderately Hard materials <i>Excellent edge retention</i>	Blue	M131	
Diallyl Phthalate (Mineral Filled) Same as Glass filled Diallyl Phthalate but, <i>More Abrasion resistant than Glass filled</i>		M130	
Transparent Powder Suited to applications that require complete visibility of the sample <i>Fair Edge retention</i>	White	M125	
Technotherm 3000 Applied across analysis that requires the a conductive sample , for SEM analysis <i>Moderate Edge retention</i>	Brown	66003630	




COLD MOUNTING SYSTEMS

PRODUCT NAME & DESCRIPTION	COMPONENTS	CODE NO
COLD MOUNTING COMPOUND Resin : 2 parts Hardener : 1 part Curing time : 15 minutes Color : Pink, opaque Used for general day to day analysis.	400grams Powder + 400 ml Liquid (1kg Weight)	MC001  

COLD MOUNTING SYSTEMS

GERMAN GRADE

PRODUCT	Mixing Ratio / Volume Ratio	Curing Time	COMPONENTS	CODE NO
Technovit Epox +(Choose A or B) <i>Transparent High clear Embedding,</i> No vacuum unit required	Resin:Hardener 2:1 		Resin (1000mL)	TECH001
		10 hours	Hardener Regular A (1000mL)	TECH002
		18 hours (at 25°C)	Hardener Fast B (1000mL)	TECH003
		3 hours (at 40°C)		
Technovit 4000 <i>Gap free, Non shrinkage embedding, white colored opaque mold</i>	Powder:Syrup1:Syrup2 2:2:1 	10 – 13 minutes	Powder (1500g) Syrup1(1000mL) Syrup 2 (500mL)	TECH026
Technovit 4002 <i>Gap free embedding, white colored opaque mold</i> Good grinding properties	Powder:Liquid 5:3 	12 – 17 minutes	Powder (1000g) Liquid (1000mL)	TECH004
Technovit 4004 <i>Transparent embedding for Routine examination</i> Vacuum unit required	Powder:Liquid 2:1 	9 – 12 minutes	Powder (1000g) Liquid (1000mL)	TECH0015

PRODUCT	Mixing Ratio / Volume Ratio	Curing Time	COMPONENTS	CODE NO
Technovit 4006 <i>High clear embedding</i> Vacuum unit required	Powder:Liquid 2:1 	9 – 13 minutes	Powder (1000g) Liquid (1000mL)	TECH0014 TECH0013
Technovit 4006 SE <i>High clear, semitransparent mold, Improved grinding properties</i> Vacuum unit required	Powder:Liquid 2:1 	11 – 15 minutes	Powder (1000g) Liquid (1000mL)	TECH030
Technovit 4071 <i>Green semitransparent mold, Quick embedding</i>	Powder:Liquid 2:1 	5 – 7 minutes	Powder (1000g) Liquid (1000mL)	TECH005
Technovit 5071 Technovit Universal Liquid <i>Green semitransparent embedding</i> <i>Dissolves in acetone</i>	Powder:Liquid 2:1 	8 – 10 minutes	Powder (1000g) Liquid (1000mL)	TECH008 TECH009
Technovit 5000 <i>Conductive material embedding, Brown opaque mold</i> <i>For SEM analysis</i>	Powder:Liquid 20g to 30ml 	7 – 12 minutes	Powder (1000g) Liquid (1000mL)	TECH023

VACUUM IMPREGNATION SYSTEM

Vacuum impregnation system is a compact pressure unit specially designed for processing transparent resins. The device is highly recommended for analysis that requires high quality bubble – free test samples.



PRODUCT NAME	Qty	CODE NO
VACUUM IMPREGNATION SYSTEM	1	VM-2000

PRECISION IMPRESSION GERMAN GRADE

For indirect inspection of the surfaces.

For 1 μ dimensional accuracy choose Technovit 3040, For <0.1 μ dimensional accuracy choose Provil Novo

PRODUCT NAME	MIXING RATIO	COMPONENTS	Qty	CODE NO
TECHNOVIT 3040	Powder:liquid ratios between 1:1 and 3:1	Technovit Universal Liquid (Choose C or D)	1000ml	TECH0010
		Yellow Powder(C)	1000g	TECH0011
		Black Powder(D)	1000g	TECH0012
PROVIL NOVO PUTTY REGULAR	Base : Catalyst 1 : 1	Base	450ml	TECH0017
		Catalyst	450ml	
PROVIL NOVO PUTTY REGULAR SOFT	Base : Catalyst 1 : 1	Base	450ml	TECH024
		Catalyst	450ml	
PROVIL NOVO LIGHT REGULAR	-	One component	2 x 50ml	TECH025



COLD MOUNTING ACCESSORIES

	PRODUCT NAME	SIZE	CODE NO
	Rubber Mould	20mm	MDSI001
		25mm	MDSI002
		30mm	MDSI003
		40mm	MDSI004
		50mm	MDSI005
		(55x25)mm	MDSI006
		(70x40)mm	MDSI007
		(100x50)mm	MDSI008
	Metallic Die	1"	MDI001
		1 ¼"	MDI002
		1 ½"	MDI003
		2"	MDI004
	Reusable Plastic Mould	30mm	MPSI011
	Plastic Holding Clips (100 nos)	0.25 (Kg)	MCL001

COLD MOUNTING ACCESSORIES

PREMIUM QUALITY *for longer life*

RUBBER MOULD:

GRADE	SIZE	CODE NO
US grade	1"	M157
	1 1/4"	M158
	1 1/2"	M159
	2"	RM2
	3"	RM3A
German grade	15mm	TECH0018
	25mm	
	30mm	
	40mm	
	50mm	

MIXING CUP AND SPATULA

PRODUCT NAME	Qty	CODE NO
Mixing Cup	10nos	TECH0020
Mixing Spatula	10nos	TECH0021

PLASTIC MOLDING CLIPS:

PRODUCT NAME	Qty/pack	CODE NO
Metal Clips	100nos	M163
Black Plastic Clips	100nos	M174




GRINDING

Grinding operation is intended to achieve flatness for further sample preparation.

The fascinating range of grinding consumables offered by Chennai Metco aids in achieving quick results



GRINDING BELTS FOR PLANAR GRINDING

PRODUCT NAME		DIMENSIONS (mm)	GRIT Sizes	CODE NO
	ABRASIVE BELT	100x915	60	BA001
	Coarse grinding of general ferrous and non ferrous materials.		80	BA002
			120	BA003
			220	BA007
			320	BA006
			400	BA004
			600	BA005
	WET BELT	100x915	80	WBA002
	Coarse grinding of general ferrous and non ferrous materials.		120	WBA003
			Can be used with water	
	ZIRCON BELT	100x915	60	ZBA001
	Coarse grinding of Stainless Steel and aluminium alloys.			
	HEAVY DUTY BELT	150X1220	60	HB001
	Coarse grinding of general ferrous and non ferrous materials.		80	HB002
			120	HB003
			Used for grinding large components.	
	LINISHING BELT	50x2000	60	LBA001
	Same as Heavy Duty belts but with different dimension		80	LBA002

SiC PAPERS FOR MANUAL GRINDING

Available Grit Size	Sheet Dimension-9"x11"
	Without PSA
80	ES001
120	ES002
180	ES003
220	ES004
320	ES005
400	ES006
600	ES007
800	ES008
1000	ES009
1200	ES010
1500	ES011
2000	ES012

SiC PAPERS FOR MANUAL GRINDING USA GRADE




Available Grit Size	Sheet Dimension 9"x11"
	Without PSA
80	M109-60
120	M109-120
180	M109-180
240	M109-240
320	M109-320
400	M109-400
600	M109-600


SiC PAPERS – PLAIN BACKING

GRINDING DISC	GRIT	8" Code no	10" Code no	12" Code no
Without PSA (Water Proof) 	80	ED001	EDE001	ESE101
	120	ED002	EDE002	ESE102
	180	ED003	EDE003	ESE103
	220	ED004	EDE004	ESE104
	320	ED005	EDE005	ESE105
	400	ED006	EDE006	ESE106
	600	ED007	EDE007	ESE107
	800	ED008	EDE008	ESE108
	1000	ED009	EDE009	ESE109
	1200	ED010	EDE010	ESE110
	1500	ED011	EDE011	ESE111
	2000	ED012	EDE012	ESE112


SiC PAPERS – PLAIN BACKING USA GRADE

GRINDING DISC	GRIT	8" Code no	10" Code no	12" Code no
Without PSA (Water Proof) 	60	M106-60	M180-60	M120-60
	120	M106-120	M180-120	M120-120
	180	M106-180	M180-180	M120-180
	240	M106-240	M180-240	M120-240
	320	M106-320	M180-320	M120-320
	400	M106-400	M180-400	M120-400
	600	M106-600	M180-600	M120-600
	800	M106-800	M180-800	M120-800
	1200	M106-1200	M180-1200	M120-1200

SiC PAPERS – PSA BACKING

GRINDING DISC	GRIT	8" Code no	10" Code no	12" Code no
With PSA (Water Proof) 	80	EA001	EAE001	ESE121
	120	EA002	EAE002	ESE122
	180	EA003	EAE003	ESE123
	220	EA004	EAE004	ESE124
	320	EA005	EAE005	ESE125
	400	EA006	EAE006	ESE126
	600	EA007	EAE007	ESE127
	800	EA008	EAE008	ESE128
	1000	EA009	EAE009	ESE129
	1200	EA010	EAE010	ESE130
	1500	EA011	EAE011	ESE131
	2000	EA012	EAE012	ESE132

SiC PAPERS – PSA BACKING USA GRADE

GRINDING DISC	GRIT	8" Code no	10" Code no	12" Code no
With PSA (Water Proof) 	60	M112-60	M160-60	M115-60
	120	M112-120	M160-120	M115-120
	180	M112-180	M160-180	M115-180
	240	M112-240	M160-240	M115-240
	320	M112-320	M160-320	M115-320
	400	M112-400	M160-400	M115-400
	600	M112-600	M160-600	M115-600
	800	M112-800	M160-800	M115-800
	1200	M112-1200	M160-1200	M115-1200

CERAMIC GRINDING FOILS – PSA

GERMAN GRADE

Suited to grinding of **Aluminium and Stainless steel alloys**

GRIT	8" Code No	10" Code No	12" Code No
80	ZA20080PSA	ZA25080PSA	ZA30080PSA
120	ZA200120PSA	ZA250120PSA	ZA300120PSA
180	ZA200180PSA	ZA250180PSA	ZA300180PSA
220	ZA200220PSA	ZA250220PSA	ZA300220PSA
320	ZA200320PSA	ZA250320PSA	ZA300320PSA
500	ZA200500PSA	ZA250500PSA	ZA300500PSA
800	ZA200800PSA	ZA250800PSA	ZA300800PSA

CERAMIC GRINDING FOILS – PET

BACKING – GERMAN GRADE

GRIT	8" Code No	10" Code No	12" Code No
80	ZA20080PET	ZA25080PET	ZA30080PET
120	ZA200120PET	ZA250120PET	ZA300120PET
180	ZA200180PET	ZA250180PET	ZA300180PET
220	ZA200220PET	ZA250220PET	ZA300220PET
320	ZA200320PET	ZA250320PET	ZA300320PET
500	ZA200500PET	ZA250500PET	ZA300500PET
800	ZA200800PET	ZA250800PET	ZA300800PET

To be used with Fast-Fix adapter discs

PREMIUM QUALITY MAGNETIC SYSTEM

Reduces the time involved in replacing PSA backed consumables

This magnetic system consists of magnetic substrate disc and a ferromagnetic carrier sheet.

PRODUCT NAME	8" Code No	10" Code No	12" Code No
MAGNETIC SUBSTRATE DISC PSA backing Pasted on the Aluminium platen	TMGD801	TMGD802	TMGD803
Ferromagnetic carrier sheet for sand paper/cloths Attached on the magnetic disc Serves as backing for PSA backed consumables	TMGDF801	TMGDF802	TMGDF803



PREMIUM QUALITY ADHESIVE SYSTEMS



PRODUCT NAME	8" Code No	10" Code No	12" Code No	14" Code No
DOUBLE PSA FOIL Double sided adhesive interface discs which are used to attach non-adhesive consumables	AS200	AS250	AS300	AS350
MD FUGO DISCS Magnetic Rubber discs with an adhesive side to attach non-adhesive papers <i>Fixation on Magnetic Platen</i>	M194-8	M194-10	M194-12	-
FAST-FIX ADAPTERDISCS SINGLE High friction surface on one side to hold the PET backed consumables. Adhesive interface on the other side to stick on to the platen.	AF200	AF250	AF300	AF350

DIAMOND EMBEDDED DISCS


Replaces the conventional silicon carbide abrasive discs.

Can process 6000 to 8000 samples. Minimizes Edge Rounding of the sample

RESIN BONDED DISCS: GERMAN GRADE


EMBEDDED DISC	SIZE μm	DIAMETER mm	CODE NO. PSA backed	CODE NO. Metal backed
APOLLO S Hard Ferrous materials 	125	200	DRS1251	DRM1251
		250	DRS1252	DRM1252
		300	DRS1253	DRM1253
		350	DRS1254	DRM1254
	91	200	DRS911	DRM911
		250	DRS912	DRM912
		300	DRS913	DRM913
		350	DRS914	DRM914
	75	200	DRS751	DRM751
		250	DRS752	DRM752
		300	DRS753	DRM753
		350	DRS754	DRM754
	54	200	DRS541	DRM541
		250	DRS542	DRM542
		300	DRS543	DRM543
		350	DRS544	DRM544
	18	200	DRS181	DRM181
		250	DRS182	DRM182
		300	DRS183	DRM183
		350	DRS184	DRM184
CARAT Fine grinding application. Nitride and oxide layer on steel can be ground without chipping. 	15	200	DRS151	DRM151
		250	DRS152	DRM152
		300	DRS153	DRM153
	6	200	DRS61	DRM61
		250	DRS62	DRM62
		300	DRS63	DRM63
	3	200	DRS31	DRM31
		250	DRS32	DRM32
		300	DRS33	DRM33

EMBEDDED DISC	SIZE μm	DIAMETER mm	CODE NO. PSA backed	CODE NO. Metal backed
SQUADRO M Patterned Surface Reduces the grinding time compared to Apollo and Carat discs	60	200	DRS601T	DRM601T
		250	DRS602T	DRM602T
		300	DRS603T	DRM603T
		350	DRS604T	DRM604T
	30	200	DRS301T	DRM301T
		250	DRS302T	DRM302T
		300	DRS303T	DRM303T
		350	DRS304T	DRM304T
	15	200	DRS151T	DRM151T
		250	DRS152T	DRM152T
		300	DRS153T	DRM153T
		350	DRS154T	DRM154T
	6	200	DRS61T	DRM61T
		250	DRS62T	DRM62T
		300	DRS63T	DRM63T
		350	DRS64T	DRM64T
	3	200	DRS31T	DRM31T
		250	DRS32T	DRM32T
		300	DRS33T	DRM33T
		350	DRS34T	DRM34T




Note :

Lubricant : Water/Water free lubricants for corrosive materials.



METAL BONDED DISCS – GERMAN GRADE

EMBEDDED DISC	SIZE μm	DIAMETER mm	CODE NO. PSA backed	CODE NO. Metal backed
PLATO High hardness ferrous material, concrete and ceramics 	125	200	DMS1251	DMM1251
		250	DMS1252	DMM1252
		300	DMS1253	DMM1253
	75	200	DMS751	DMM751
		250	DMS752	DMM752
		300	DMS753	DMM753
	54	200	DMS541	DMM541
		250	DMS542	DMM542
		300	DMS543	DMM543
	10	200	DMS101	DMM101
		250	DMS102	DMM102
		300	DMS103	DMM103



Note :

Lubricant : Water/Water free lubricants for corrosive materials.

RESIN BONDED DISCS: USA GRADE

Flexible Magnetic Backing



GRIT	8" (203mm)	10"(254mm)	12"(305mm)
80	CM506-8G	CM506-10G	CM506-12G
120	CM506-8R	CM506-10R	CM506-12R
220	CM506-8Y	CM506-10Y	CM506-12Y
600	CM506-8B	CM506-10B	CM506-12B
1200	CM506-8O	CM506-10O	CM506-12O



METAL BONDED DISCS: USA GRADE

PSA backing

Micron	8" (203mm)	12"(305mm)
165	M190-165	M191-165
70	M190-70	M191-70
45	M190-45	M191-45
30	M190-30	M191-30
15	M190-15	M191-15
9	M190-9	M191-9



SPECTRO GRINDING DISCS

Huge samples subjected to spectrographic analysis are coarse ground through these discs. The disc is composed of coarse alumina abrasive.

GRIT SIZE	CODE NO	SIZE (Dia x ID) mm	BACKING TYPE
36	SGD002	350X40	Without PSA
60	SGD003		
80	SGD004		
40	SGD006PS	350X40	With PSA




PORTABLE GRINDING DISCS

These grinding discs are used in **Portable Grinding Machines**.

Applications that require the work piece to be tested in its original environment. Example: foundries, welding shops, fabrication plants, steel mills, forge shops and shipyard.

Available Grit Size (30mm)	Code No With PSA
80	SED007
120	SED008
180	SED009
220	SED010
320	SED011
400	SED012
600	SED013
800	SED014
1000	SED015
1200	SED016
1500	SED017
2000	SED018





POLISHING

Polishing operation is aimed at achieving mirror finish on the sample surface for further microscopic analysis. The right combination of material with proper polishing media yields a perfect finish most suitable for the application.

Chennai Metco offers a vast range of polishing consumables for most demanding applications.

POLISHING CLOTHS

CLOTH	NAP	FABRIC	RECOMMENDED USE
BILLIARDS <i>Rough Polishing</i>	Napless 	Wool	Used with 6 – 12 μ diamond paste and Aerosol spray.
VELVET <i>Fine Polishing</i>	Long napped 	Synthetic	Used with 3 – 4 μ diamond paste and Aerosol spray.
MICRO CLOTH <i>Ultra-Fine Polishing</i>	Short Nap 	Synthetic rayon	Used with 0.5 – 1 μ diamond paste and Aerosol spray.
SELVYT <i>For Perfect mirror finish</i>	Short Piled 	Cotton	Alumina / colloidal silica suspensions.





POLISHING CLOTH – WITHOUT PSA



POLISHING CLOTH	PRODUCT NAME	Qty	10" Code No	12" Code No	14" Code No
Without PSA	Velvet	10	PC001	PC12V	PC14V
	Micro	10	PC002	PC12M	PC14M
	Selvyt	10	PC003	PC12S	PC14S
	Billiards	10	PC004	PC12B	PC14B
	Velvet (Black)	10	LP10PC	LP12PC	
	Velvet (Brown)	10	LP10BR	LP12BRPC	

POLISHING CLOTH – WITH PSA

POLISHING CLOTH	PRODUCT NAME	Qty	8" Code No	10" Code No
With PSA	Velvet Black	10	LP001	LP002
	Velvet Brown	10	LP003	LP004

POLISHING CLOTHS – USA GRADE

CLOTH	Cloth Design	NAP	FABRIC	RECOMMENDED USE
CHEM PAD Final polishing on titanium, stainless, lead/tin solders, electronic packages and soft non-ferrous and plastics		Napless	Soft synthetic chemically resistant	Used with aluminum oxide powders, colloidal silica suspensions and 1 μ and finer diamond products.
DIA DAC Pre-polishing stages on ferrous and non-ferrous materials, petrographic samples and ceramics.		Napless	Satin acetate woven pure medium hard silk cloth of very close weave	Used with 3 μ diamond products.
DIA PLUS Polishing on Ferrous and sintered carbides. Excellent for unmounted samples or manual polishing		Long Napped	Synthetic Cloth	Used with 3 μ and finer diamond products, alumina suspensions.
DIA PLAN Maintain Flatness on hard / very hard materials		Napless	Woven polyester nylon cloth of close weave	Used with 15 to 3 μ diamond products.
DIA MOL Typically used for softer materials		Napless	Woven wool medium / soft	Used with 6 to 1 μ diamond. Alumina powders.

CLOTH	Cloth Design	NAP	FABRIC	RECOMMENDED USE
DIA-NAP Final Polishing of all ferrous metals and sintered carbides		very short napped	synthetic cloth	Used with 1 μ and finer diamond
VEL LAP Used for final polishing of soft metals		Short Nap	white synthetic	Used with Aluminas, diamond and colloidal silica suspensions.
MICRO CLOTH Final polishing of all materials Higher removal		Short Nap	synthetic rayon	Used with Alumina and 1 μ diamond suspension
RED FELT Rough and intermediate polishing of all metals		Napless	100% pure virgin wool fabric plucked rather than sheared pile	Used with 1 μ diamond and Alumina suspensions
PAN W Maintain flatness, edge and inclusion retention in pre-polishing stages on metals, ceramics etc		Napless	Hard Synthetic Chemotextile	Used with 3 μ or finer diamond compounds


POLISHING CLOTHS WITH PSA USA GRADE


<i>With PSA</i>	NOS / Pack	8" CODE NO	10" CODE NO	12" CODE NO
Dia-Mol	5	M620-8PS	M620-10PS	M620-12PS
Dia Nap	5	M630-8PS	M630-10PS	M630-12PS
Chem Pad	5	M660-8PS	M660-10PS	M660-12PS
Dia Dac	5	M615-8PS	M615-10PS	M615-12PS
Dia Plan	5	M610-8PS	M610-10PS	M610-12PS
Dia Plus	5	M625-8PS	M625-10PS	M625-12PS
Red Felt	10	M570-8PS	M570-10PS	M570-12PS
Micro	10	M500-8PS	M500-10PS	M500-12PS
Vel Lap	10	M536-8PS	M536-10PS	M536-12PS
Pan w	10	M560-8PS	M560-10PS	M560-12PS

POLISHING CLOTHS MAGNETIC BACKING - USA GRADE

<i>Magnetic Backing</i>	NOS / Pack	8" CODE NO	10" CODE NO	12" CODE NO
Dia-Mol	5	M620-8MD	M620-10MD	M620-12MD
Pan w	5	M560-8MD	M560-10MD	M560-12MD
Micro	5	M500-8MD	M500-10MD	M500-12MD
Dia Nap	5	M630-8MD	M630-10MD	M630-12MD
Chem Pad	5	M660-8MD	M660-10MD	M660-12MD
Dia Dac	5	M615-8MD	M615-10MD	M615-12MD
Dia Plan	5	M610-8MD	M610-10MD	M610-12MD
Dia Plus	5	M625-8MD	M625-10MD	M625-12MD

POLISHING CLOTHS – GERMAN GRADE

CLOTH		THICKNESS	FABRIC	RECOMMENDED USE
QUICK-STEP First choice for pre-polishing . Very high removal rate on all materials		0.5 mm	Hard synthetic fibre	Used with 6 - 15µ diamond products
STEP-PLUS Good removal rate and superb finishing on almost all materials		0.75 mm	Medium hard synthetic fiber cloth	Used with 6 - 3µ diamond products
SAMBA-N Fine polishing of larger workpieces, excellent for manual polishing		0.65 mm	Medium hard woven wool polishing cloth	Used with 1 - 15µ diamond products
JAZZ-PLUS Good removal rate and superb finish on metals and optical materials preparation		0.7 mm	Non-woven polishing pad, PU impregnated	Used with 1 - 6µ diamond products

CLOTH		THICKNESS	FABRIC	RECOMMENDED USE
FOX-PLUS Good removal rate and superb finishing on all kinds of metals		0.85 mm	Flocked universal polishing pad	Used with 1 - 15μ diamond products
SWING-PLUS Polishing pad for super-finishing.		0.7 mm	Flocked universal polishing cloth with very low flock density	Used with 0.1 - 6μ diamond products
ALUPOL-PLUS Perfect cross-over pad, for diamond and oxide polishing.		0.6 mm	Flocked universal polishing cloth with extremely high flock density	Used with 0.1 - 3μ diamond, Alumina suspensions
MAMBO Final finishing with oxide polishing slurries		1.5 mm	Chemically resistant, spongy oxide polishing cloth	Alumina and colloidal silica suspensions

POLISHING CLOTHS WITH PSA GERMAN GRADE

<i>With PSA</i>	8" CODE NO	10" CODE NO	12" CODE NO	14" CODE NO	16" CODE NO
Quick-Step	PCG01	PCG11	PCG21	PCG31	PCG41
Step-Plus	PCG02	PCG12	PCG22	PCG32	PCG42
Samba-N	PCG03	PCG13	PCG23	PCG33	PCG43
Jazz-Plus	PCG04	PCG14	PCG24	PCG34	PCG44
Fox-Plus	PCG05	PCG15	PCG25	PCG35	PCG45
Swing-Plus	PCG06	PCG16	PCG26	PCG36	PCG46
Alupol-Plus	PCG07	PCG17	PCG27	PCG37	PCG47
Mambo	PCG08	PCG18	PCG28	PCG38	PCG48

POLISHING CLOTHS METAL BACKING GERMAN GRADE

<i>Metal Backing</i>	8" CODE NO	10" CODE NO	12" CODE NO	14" CODE NO
Quick-Step	PDG01	PDG11	PDG21	PDG31
Step-Plus	PDG02	PDG12	PDG22	PDG32
Samba-N	PDG03	PDG13	PDG23	PDG33
Jazz-Plus	PDG04	PDG14	PDG24	PDG34
Fox-Plus	PDG05	PDG15	PDG25	PDG35
Swing-Plus	PDG06	PDG16	PDG26	PDG36
Alupol-Plus	PDG07	PDG17	PDG27	PDG37
Mambo	PDG08	PDG18	PDG28	PDG38

DIAMOND SLURRIES - USA GRADE

Polycrystalline *Water based* diamond slurries

Should be used with Lubricants

μm	Bottle only 8 oz (240ml)	Air Spray 8 oz (240ml)	Trigger Spray 16 oz (475ml)	Refill bottle 32 oz (950ml)	Refill bottle 128 oz (3.8L)
0.25	M270	M280	M280-16	M280-32	M280-128
1	M271	M281	M281-16	M281-32	M281-128
3	M272	M282	M282-16	M282-32	M282-128
6	M273	M283	M283-16	M283-32	M283-128
9	M274	M284	M284-16	M284-32	M284-128
15	M275	M285	M285-16	M285-32	M285-128
30	M276	M286	M286-16	M286-32	M286-128
45	M277	M287	M287-16	M287-32	M287-128

Monocrystalline *Water based* diamond slurries

Should be used with Lubricants

μm	Bottle only 8 oz (240ml)	Finger pump 8 oz (240ml)	Trigger Spray 16 oz (475ml)	Refill bottle 32 oz (950ml)	Refill bottle 128 oz (3.8L)
0.10	M250	M240	M260	M260-32	M260-128
0.25	M251	M241	M261	M261-32	M261-128
0.5	M252	M242	M262	M262-32	M262-128
1	M253	M243	M263	M263-32	M263-128
3	M254	M244	M264	M264-32	M264-128
6	M255	M245	M265	M265-32	M265-128
9	M256	M246	M266	M266-32	M266-128
15	M257	M247	M267	M267-32	M267-128
30	M258	M248	M268	M268-32	M268-128
45	M259	M249	M269	M269-32	M269-128



Polycrystalline *Oil based* diamond slurries

Should be used with Lubricants.

Suited to water sensitive materials

Micron	Bottle only 8 oz (240ml)	Trigger Spray 16 oz (475ml)
0.25	M260-P25	M260-16
1	M261-P1	M261-16
3	M262-P3	M262-16
6	M260-P6	M260-P-16
9	M261-P9	M261-P-16
15	M262-P15	M262-P-16

Lubricants to be used with diamond slurries
(USA Grade)

Type	Quart (950ml)	Gallon (3.8L)	5 Gallon (19L)
Purple (Alcohol based)	M225-32	M225-128	M225-640
Green (Water based)	M226-32	M226-128	M226-640
Red (Water/Oil based)	M227-32	M227-128	M227-640
Brown (Alcohol based)	M220-32	-	-
Blue Lube (Ethanol based)	M222	M223	M224



Lubricant Integrated Diamond Slurries

USA GRADE

Description	Compatible Cloth	16 oz (475ml)	128 oz (3.8L)
DiaLube Plan	Dia Plan	M281-16DL	M281-128DL
DiaLube Mol	Dia Mol	M282-16DL	M282-128DL
DiaLube Dac	Dia Dac	M283-16DL	M283-128DL
DiaLube Plus	Dia Plus	M284-16DL	M284-128DL
DiaLube Nap	Dia Nap	M285-16DL	M285-128DL

A unique blend of quality high performance diamonds accompanied with cooling lubricants specially formulated to suit the respective polishing cloths. Can be used without lubricants

DIAMOND SLURRIES - GERMAN GRADE



MAGNUM-TOP **Monocrystalline** *Water based* diamond slurries

µm	250ml	500ml	1L	2.5L	10L
0.1	DS100W	DS200W	DS300W	DS400W	DS500W
0.25	DS101W	DS201W	DS301W	DS401W	DS501W
1	DS102W	DS202W	DS302W	DS402W	DS502W
3	DS103W	DS203W	DS303W	DS403W	DS503W
6	DS104W	DS204W	DS304W	DS404W	DS504W
9	DS105W	DS205W	DS305W	DS405W	DS505W
15	DS106W	DS206W	DS306W	DS406W	DS506W
30	DS107W	DS208W	DS307W	DS407W	DS507W
45	DS108W	DS208W	DS308W	DS408W	DS508W

Oil based suspensions for water sensitive materials available

POLY-TOP **Polycrystalline** *Water based* diamond slurries

µm	250ml	500ml	1L	2.5L	10L
0.1	DS150W	DS250W	DS350W	DS450W	DS550W
0.25	DS151W	DS251W	DS351W	DS451W	DS551W
1	DS152W	DS252W	DS352W	DS452W	DS552W
3	DS153W	DS253W	DS353W	DS453W	DS553W
6	DS154W	DS254W	DS354W	DS454W	DS554W
9	DS155W	DS255W	DS355W	DS455W	DS555W
15	DS156W	DS256W	DS356W	DS456W	DS556W

Oil based suspensions for water sensitive materials available

Lubricants to be used with diamond slurries
(German Grade)

Type	1L	5L	10L	200L
Lubricant Yellow (Water base)	LY1	LY2	LY3	-
Lubricant Yellow concentrate	-	-	LYC3	-
Lubricant Blue (Alcohol base)	LB1	LB2	LB3	LB4
Lubricant Blue concentrate	LBC1	-	-	-
Lubricant Orange (Oil base)	LO1	-	LO3	LO4
Lubricant Red (Water/oil base)	LR1	-	LR3	-
Lubricant White (Alcohol base)	LW1	-	LW3	-
Lubricant White concentrate	LWC1	-	-	-

Water base diamond Slurries + Lubricants ↓ Quick Results with just a few drops of slurry and Lubricant	Lubricant Integrated Slurries ↓ Results achieved at a moderate rate
---	--

Note: The above table is applicable for German grade diamond slurries

Lubricant Integrated Diamond Slurries

GERMAN GRADE

MAGNUM-TOP DUO **Monocrystalline**

Lubricant Integrated **diamond slurries**

Can be used without lubricants

µm	250ml	500ml	1L	2.5L	10L
0.1	DS100L	DS200L	DS300L	DS400L	DS500L
0.25	DS101L	DS201L	DS301L	DS401L	DS501L
1	DS102L	DS202L	DS302L	DS402L	DS502L
3	DS103L	DS203L	DS303L	DS403L	DS503L
6	DS104L	DS204L	DS304L	DS404L	DS504L
9	DS105L	DS205L	DS305L	DS405L	DS505L
15	DS106L	DS206L	DS306L	DS406L	DS506L
20	DS107L	DS208L	DS307L	DS407L	DS507L
30	DS108L	DS208L	DS308L	DS408L	DS508L

POLY-TOP DUO **Polycrystalline** *Lubricant*

Integrated **diamond slurries**

Can be used without lubricants

µm	250ml	500ml	1L	2.5L	10L
0.1	DS150L	DS250L	DS350L	DS450L	DS550L
0.25	DS151L	DS251L	DS351L	DS451L	DS551L
1	DS152L	DS252L	DS352L	DS452L	DS552L
3	DS153L	DS253L	DS353L	DS453L	DS553L
6	DS154L	DS254L	DS354L	DS454L	DS554L
9	DS155L	DS255L	DS355L	DS455L	DS555L
15	DS156L	DS256L	DS356L	DS456L	DS556L

DIAMOND PASTE

	Product	Particle Size	Weight	CODE NO
	Diamond Paste	0.5-1 Micron	0.005Kg/5cc	DP001
		3-4 Micron		DP002
		6-12 Micron		DP003

LUBRICANT AEROSOL SPRAY

	PRODUCT	QTY	CODE NO
	Lubricant Aerosol Spray	200ml	DL001

DIAMOND PASTE – USA GRADE



Product	μm	CODE NO 5g	CODE NO 10g	CODE NO 20g
MONOCRYSTALLINE Water and Oil Soluble	0.25	M200	-	M201
	1	M202	-	M203
	3	M204	-	M205
	6	M206	-	M207
	9	M208	-	M209
	15	M210	-	M211
	30	M212	-	M213
	45	M214	-	M215
POLYCRYSTALLINE Oil Soluble	0.25	-	M200P	-
	1	-	M202P	-
	3	-	M204P	-
	6	-	M206P	-
	9	-	M208P	-
	15	-	M210P	-

DIAMOND AEROSOL SPRAY USA GRADE

PRODUCT	μm	QTY	CODE NO
DIAMOND AEROSOL SPRAY 	0.25	150ml	M288
	1		M289
	3		M290
	6		M291
	9		M292
	15		M293
	30		M294
	45		M295

DIAMOND PASTE – GERMAN GRADE



Diamond Paste – Universal (Oil, Water and Alcohol soluble)

No Lubricant needed

Product	µm	Code No 5g	Code No 10g	Code No 20g
SPARKLING U/M	0.25	DP100U	DP200U	DP300U
Monocrystalline	1	DP101U	DP201U	DP301U
	3	DP102U	DP202U	DP302U
	6	DP103U	DP203U	DP303U
	9	DP104U	DP204U	DP304U
	15	DP105U	DP205U	DP305U
	20	DP106U	DP206U	DP306U
	30	DP107U	DP207U	DP307U
	45	DP108U	DP208U	DP308U
	70	DP109U	DP209U	DP309U
	90	DP110U	DP210U	DP310U
SPARKLING U/P	0.25	DP150U	DP250U	DP350U
Polycrystalline	1	DP151U	DP251U	DP351U
	3	DP152U	DP252U	DP352U
	6	DP153U	DP253U	DP353U
	9	DP154U	DP254U	DP354U
	15	DP155U	DP255U	DP355U
	20	DP156U	DP256U	DP356U
	30	DP157U	DP257U	DP357U

Diamond Paste – Oil Soluble

**Corrosion-Free polishing*

**No Lubricant needed*

**Both Automatic and Manual polishing*

Product	µm	Code No 5g	Code No 10g	Code No 20g
SPARKLING O/M Monocrystalline	0.25	DP1000	DP2000	DP3000
	1	DP1010	DP2010	DP3010
	3	DP1020	DP2020	DP3020
	6	DP1030	DP2030	DP3030
	9	DP1040	DP2040	DP3040
	15	DP1050	DP2050	DP3050
	20	DP1060	DP2060	DP3060
	30	DP1070	DP2070	DP3070
	45	DP1080	DP2080	DP3080
	70	DP1090	DP2090	DP3090
	90	DP1100	DP2100	DP3100
SPARKLING O/P Polycrystalline	0.25	DP1500	DP2500	DP3500
	1	DP1510	DP2510	DP3510
	3	DP1520	DP2520	DP3520
	6	DP1530	DP2530	DP3530
	9	DP1540	DP2540	DP3540
	15	DP1550	DP2550	DP3550
	20	DP1560	DP2560	DP3560
	30	DP1570	DP2570	DP3570

Diamond Paste – Oil and Alcohol Soluble

**Preferred for steels*

**No Lubricant needed*

**Both Automatic and Manual polishing*

Product	µm	Code No 5g	Code No 10g	Code No 20g
SPARKLING FAS/M Monocrystalline	0.25	DP1000A	DP2000A	DP3000A
	1	DP1010A	DP2010A	DP3010A
	3	DP1020A	DP2020A	DP3020A
	6	DP1030A	DP2030A	DP3030A
	9	DP1040A	DP2040A	DP3040A
	15	DP1050A	DP2050A	DP3050A
	20	DP1060A	DP2060A	DP3060A
	30	DP1070A	DP2070A	DP3070A
	45	DP1080A	DP2080A	DP3080A
	70	DP1090A	DP2090A	DP3090A
	90	DP1100A	DP2100A	DP3100A
SPARKLING FAS/P Polycrystalline	0.25	DP1500A	DP2500A	DP3500A
	1	DP1510A	DP2510A	DP3510A
	3	DP1520A	DP2520A	DP3520A
	6	DP1530A	DP2530A	DP3530A
	9	DP1540A	DP2540A	DP3540A
	15	DP1550A	DP2550A	DP3550A
	20	DP1560A	DP2560A	DP3560A
	30	DP1570A	DP2570A	DP3570A

Diamond Paste – Water Soluble

**Easy Cleaning*

**Used with Lubricant Yellow*

**Both Automatic and Manual polishing*

Product	µm	Code No 5g	Code No 10g	Code No 20g
SPARKLING W/M	0.25	DP100W	DP200W	DP300W
Monocrystalline	1	DP101W	DP201W	DP301W
	3	DP102W	DP202W	DP302W
	6	DP103W	DP203W	DP303W
	9	DP104W	DP204W	DP304W
	15	DP105W	DP205W	DP305W
	20	DP106W	DP206W	DP306W
	30	DP107W	DP207W	DP307W
	45	DP108W	DP208W	DP308W
	70	DP109W	DP209W	DP309W
	90	DP110W	DP210W	DP310W
SPARKLING W/P	0.25	DP150W	DP250W	DP350W
Polycrystalline	1	DP151W	DP251W	DP351W
	3	DP152W	DP252W	DP352W
	6	DP153W	DP253W	DP353W
	9	DP154W	DP254W	DP354W
	15	DP155W	DP255W	DP355W
	20	DP156W	DP256W	DP356W
	30	DP157W	DP257W	DP357W

***Note: Diamond Paste in Stick form available
Bulk quantity (400g and 1000g) available in
Jar form***

DIAMOND AEROSOL SPRAY GERMAN GRADE

Product	Micron	200ml
POLY DIAMOND SPRAY Polycrystalline	0.25	DA150W
	1	DA151W
	3	DA152W
	6	DA153W
	9	DA154W
	15	DA155W

ALUMINA COMPOUNDS

For Ultra-Fine polishing

Product	µm	Qty	CODE NO
Alumina Powder	1	1 Kg	AP001
Alumina Fast Pol (Alumina Suspension)	1	200ml	IFP00001



PREMIUM QUALITY ALUMINA AND COLLOIDAL SILICA COMPOUNDS

For Ultra-Fine polishing

De-agglomerated alumina



Product	μm	Qty	CODE NO
Alumina Powder	0.05	455g	M300
	0.3	455g	M301
	1	455g	M302
Alumina Suspensions	0.05	200ml	IFP00004
	0.1	1000ml	M310-32
	0.3	200ml	IFP00003
	1	200ml	IFP00001
Colloidal Silica <i>pH 10.5</i>	0.04	1000ml	M511-OPSQ
	0.05	1L	OPS00001
	0.06	1000ml	M306-32


PORTABLE POLISHING CLOTH

The portable polishing cloths find great application in **portable polishing devices**. Applications that involve **on-site inspection of metals** require these portable cloths to polish the surface of interest.

POLISHING CLOTH	PRODUCT NAME	SIZE(mm)	CODE NO
With PSA	Velvet Polishing Cloth-Grey	30	PCS002
	Velvet Polishing Cloth-Black(Premium Quality)		LP30PS
	Velvet Polishing Cloth Brown (Premium Quality)		LP30BRPS
With PSA (German Grade)	Technotron MM 420	75	TECH029

ACCESSORIES

	PRODUCT	QUANTITY	CODE NO.
	Hot Air Dryer	1	HTP001
	Specimen Leveller	1	SPL001
	Portable Engraving Machine	1	EG001
	Vacuum Impregnation System	1	VM-2000
	Specimen Leveling Clay	200gm	SLC001
	Sample Protective Liquid (Dissolves in acetone)	200ml	LA001



Chennai Metco Pvt Ltd.,

SP100A, Ambattur Industrial Estate, Chennai - 600 058.

Ph: 91-044-2688 0005, 91-044-3999 5000, fax: 91-044-2688 0850

e-mail: mktg@chennai metco.com, mktg@chennai metco.co.in

www.chennai metco.com



This document was created with Win2PDF available at <http://www.win2pdf.com>.
The unregistered version of Win2PDF is for evaluation or non-commercial use only.
This page will not be added after purchasing Win2PDF.